
Product Data Sheet

Product Name: Alexamorelin Met 1

Cat. No.: GC35282

Chemical Properties

Cas. No.

Formula

M.Wt

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Alexamorelin Met 1 is one of the metabolites of alexamorelin. The heptapeptide Ala-His-D-2-methyl-Trp-Ala-Trp-D-Phe-Lys-NH₂ (Alexamorelin) is a synthetic molecule which inhibits growth hormone secretagogue binding in vitro.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA