
Product Data Sheet

Product Name: Alosetron D3 Hydrochloride

Cat. No.: GC35303

Chemical Properties

Cas. No. 1189919-71-8

SMILES O=C1N(CC2=C(C)N=CN2)CCC3=C1C4=CC=CC=C4N3C([2H])([2H])[2H].[H]Cl

Formula $C_{17}H_{16}D_3ClN_4O$ M.Wt 333.83

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Alosetron D3 Hydrochloride (GR-68755C D3) is deuterium labeled Alosetron, which is a serotonin 5HT₃-receptor antagonist.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA