

---

**Product Data Sheet**

---

Product Name: Angelic anhydride

Cat. No.: GC35344

**Chemical Properties**

Cas. No. 94487-74-8

SMILES C/C=C(C)\C(OC(/C(C)=C\C)=O)=OFormula  $C_{10}H_{14}O_3$  M.Wt 182.22

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Angelic anhydride (2-Methylisocrotonic Anhydride, (Z)-2-Methylbut-2-enoic anhydride) is a synthetic chemical and usually used in content determination assays and pharmacological experiments.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA