
Product Data Sheet

Product Name: C-Reactive Protein (CRP) 174-185

Cat. No.: GC35744

Chemical Properties

Cas. No. 160369-86-8

Formula $C_{62}H_{93}N_{13}O_{16}$ M.Wt 1276.48

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

C-Reactive Protein (CRP) 174-185 is the 174-185 fragment of C-Reactive Protein. C-Reactive Protein (CRP), the prototypic marker of inflammation, is a cardiovascular risk marker and may promote atherogenesis[1].

[1]. Devaraj S, et al. C-reactive protein stimulates superoxide anion release and tissue factor activity in vivo. *Atherosclerosis*. 2009 Mar;203(1):67-74.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
