
Product Data Sheet

Product Name: Dye 937
 Cat. No.: GC35915

Chemical Properties

Cas. No. 195199-04-3

SMILES C[N+]1=C(/C=C2C(C=CC=C3)=C3N(C4=CC=CC=C4)C(N(CCC)CCCN(C)C)=C\2)SC5=CC=CC=C51.[I-]

Formula $C_{32}H_{37}N_4S$ M.Wt 636.63

Solubility Soluble in DMSO Storage Store at -20°C

General For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice Condition upon request.

Structure

Background

Dye 937, substituted unsymmetrical cyanine dyes with selected permeability, useful in the detection of DNA in electrophoretic gels.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
