
Product Data Sheet

Product Name: E3 ligase Ligand 18

Cat. No.: GC35932

Chemical Properties

Cas. No. 2241669-88-3

SMILES N[C@H](C(C)(C)C)C(N1[C@H](C(NCC2=CC=C(C3=C(C)N=CS3)C=C2)=O)C[C@@H](O)C1)=O

Formula	C ₂₂ H ₃₀ N ₄ O ₃ S	M.Wt	430.56
---------	---	------	--------

Solubility	Soluble in DMSO	Storage	Store at -20°C
------------	-----------------	---------	----------------

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

E3 ligase Ligand 18 is a ligand for E3 ubiquitin ligase. E3 ligase Ligand 18 can be connected to the ligand for protein (e.g., BCR-ABL1) by a linker to form PROTACs (e.g., GMB-475). GMB-475 induces the degradation of BCR-ABL1 with an IC₅₀ of 1.11 μM in Ba/F3 cells[1]. VHL

[1]. Burslem GM, et al. Targeting BCR-ABL1 in Chronic Myeloid Leukemia by PROTAC-mediated Targeted Protein Degradation. Cancer Res. 2019 Jul 16. pii: canres.1236.2019.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
