
Product Data Sheet

Product Name: E3 ligase Ligand-Linker Conjugates 40

Cat. No.: GC35947

Chemical Properties

Cas. No. 2222354-20-1

SMILES C[C@H](N(C(OC(C)(C)C)=O)C)C(N[C@@H](C1CCCCC1)C(N(CCC2[C@@H]2C3=NC(C(C4=CC=CC(OCCOCCOCCOCC(O)=O)=C4)=O)=CS3)=O)=O

Formula	C ₃₉ H ₅₆ N ₄ O ₁₁ S	M.Wt	788.95
---------	--	------	--------

Solubility	Soluble in DMSO	Storage	Store at -20°C
------------	-----------------	---------	----------------

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

E3 ligase Ligand-Linker Conjugates 40 incorporates a cIAP ligand for the E3 ubiquitin ligase, and a PROTAC linker. E3 ligase Ligand-Linker Conjugates 40 can be used to design PROTAC degrader[1]. cIAP1

[1]. Itoh Y, et al. Design, synthesis and biological evaluation of nuclear receptor-degradation inducers. *Bioorg Med Chem.* 2011 Nov 15;19(22):6768-78.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
