
Product Data Sheet

Product Name: F1324
Cat. No.: GC36022

Chemical Properties

Cas. No. 2607975-44-8

Formula M.Wt

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

F1324 is a potent, high affinity peptidic inhibitor of B cell lymphoma 6 (BCL6) with an IC50 of 1 nM[1]. IC50: 1 nM (BCL6)[1]

[1]. Sakamoto K, et al. Discovery of high-affinity BCL6-binding peptide and its structure-activity relationship. Biochem Biophys Res Commun. 2017 Jan 8;482(2):310-316.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA