
Product Data Sheet

Product Name: Ganoderenic acid E

Cat. No.: GC36115

Chemical Properties

Cas. No. 110241-23-1

SMILES CC12C3=C(C(C(O)C1(C/C(C)=C\C(CC(C)C(O)=O)=O)CC2=O)C)=O)C4(C(C(C)(C(CC4)=O)C)CC3O)CFormula C₃₀H₄₀O₈ M.Wt 528.63

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Ganoderenic acid E is a lanostane-type triterpene isolated from Ganoderma lucidum.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA