
Product Data Sheet

Product Name: Glycoprotein 276-286

Cat. No.: GC36162

Chemical Properties

Cas. No. 160543-97-5

Formula C₄₆H₇₀N₁₂O₁₇S M.Wt 1095.18

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Glycoprotein (276-286) is a Db-restricted peptide derived from lymphocytic choriomeningitis virus (LCMV) glycoprotein (GP), corresponds to amino acids 276-286[1].

[1]. Oldstone MB, et al. Discriminated selection among viral peptides with the appropriate anchor residues: implications for the size of the cytotoxic T-lymphocyte repertoire and control of viral infection. J Virol. 1995 Dec;69(12):7423-9.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
