
Product Data Sheet

Product Name: L-Leucyl-L-alanine

Cat. No.: GC36412

Chemical Properties

Cas. No. 7298-84-2

Formula $C_9H_{18}N_2O_3$

M.Wt 202.25

Solubility Water: 100 mg/mL (494.44 mM); Methanol: 8.33 mg/mL
(41.19 mM)Storage Store at -
20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

L-Leucyl-L-alanine Hydrate (H-Leu-Ala-OH) is a dipeptide composed of L-leucine and L-alanine joined by a peptide linkage. It is a metabolite.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA