
Product Data Sheet

Product Name: MC-VC-PABC-DNA31

Cat. No.: GC36577

Chemical Properties

Cas. No. 1639352-03-6

Formula C₇₇H₉₄N₁₀O₂₁ M.Wt 1495.62

Solubility Soluble in DMSO Storage 4°C, away from moisture and light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

MC-VC-PABC-DNA31 has a bioreversible linkage based on a quaternary ammonium for targeted delivery and it can improve pharmacokinetics and the therapeutic index. MC-VC-PABC-DNA31 is used for the antibody-drug conjugates (ADC) that are effective and stable in vitro and in vivo to treat various diseases or disorders[1].

[1]. Staben LR, et al. Targeted drug delivery through the traceless release of tertiary and heteroaryl amines from antibody-drug conjugates. Nat Chem. 2016 Dec;8(12):1112-1119.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA