
Product Data Sheet

Product Name: Mini Gastrin I, human

Cat. No.: GC36610

Chemical Properties

Cas. No. 54405-27-5

Formula $C_{74}H_{99}N_{15}O_{26}S$ M.Wt 1646.73

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

Mini Gastrin I, human is a shorter version of human gastrin, consists of amino acids 5-17 of the parent peptide. CCK2i4svR[1]

[1]. Laverman P, et al. Targeting of a CCK(2) receptor splice variant with (111)In-labelled cholecystokinin-8 (CCK8) and (111)In-labelled minigastrin. Eur J Nucl Med Mol Imaging. 2008 Feb;35(2):386-92.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA