
Product Data Sheet

Product Name: N-acetyl Dapsone D4'

Cat. No.: GC36683

Chemical Properties

Cas. No.

SMILES #583T2OrInOmFkOrYB+VMfQudA5NfA2UFkOrkNioINrTBVoT2OoTNWbdkOoTBVUFk/UTNWeOnWSdlquOj5NfA2UFk/=

Formula C₁₄H₁₀D₄N₂O₃S

M.Wt

294.36

Solubility Soluble in DMSO

Storage

Store at -20°C

General For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock tips solution can be stored below -20°C for several months.

Shipping Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.
Condition

Structure

Background

N-acetyl Dapsone (D4') is the deuterium labeled N-acetyl Dapsone, which is a metabolite of Dapsone.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA