
Product Data Sheet

Product Name: OVA sequence 323-336

Cat. No.: GC36827

Chemical Properties

Cas. No. 92915-80-5

Formula $C_{63}H_{100}N_{20}O_{22}$

M.Wt

1489.59

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

OVA sequence (323-336) is a cognate helper T-lymphocyte peptide that is employed to enhance CTL epitope immunogenicity.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA