
Product Data Sheet

Product Name: PSA1 141-150

Cat. No.: GC37023

Chemical Properties

Cas. No. 160215-60-1

Formula C₅₅H₉₃N₁₃O₁₃S M.Wt 1176.47

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

PSA1 (141-150), a prostate specific antigen 1 peptide, is used in the immunotherapy of cancer experiments[1].

[1]. Arlen PM, et al. Pox viral vaccine approaches. Semin Oncol. 2005 Dec;32(6):549-55.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
