
Product Data Sheet

Product Name: p-Syneprine

Cat. No.: GC37033

Chemical Properties

Cas. No. 614-35-7

SMILES OC1=CC=C([C@@H](O)CNC)C=C1Formula $C_9H_{13}NO_2$

M.Wt 167.21

Solubility Water: 2 mg/mL (11.96 mM)

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

p-Syneprine is an organic compound, found in multiple biofluids, such as urine and blood. Human Endogenous Metabolite

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA