

---

**Product Data Sheet**

---

Product Name: SIYRY  
Cat. No.: GC37645

**Chemical Properties**

Cas. No. 178561-37-0

Formula C<sub>50</sub>H<sub>71</sub>N<sub>11</sub>O<sub>13</sub> M.Wt 1034.16

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

**Background**

SIYRY is a Kb-restricted epitope peptide[1].

[1]. Sykulev Y, et al. Peptide antagonism and T cell receptor interactions with peptide-MHC complexes. *Immunity*. 1998 Oct;9(4):475-83.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA