

---

**Product Data Sheet**

---

Product Name: Super-TDU (1-31) TFA

Cat. No.: GC37704

**Chemical Properties**

Cas. No.

Formula  $C_{141}H_{218}N_{40}O_{48}.C_2HF_3O_2$  M.Wt 3355.5Solubility Soluble in DMSO Storage  $-20^{\circ}C$ , protect from light

General tips For obtaining a higher solubility , please warm the tube at  $37^{\circ}C$  and shake it in the ultrasonic bath for a while. Stock solution can be stored below  $-20^{\circ}C$  for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Super-TDU (1-31) is a peptide of Super-TDU, which is an inhibitor of YAP-TEADs, shows potent anti-tumor activity[1].

[1]. Jiao S, et al. A peptide mimicking VGLL4 function acts as a YAP antagonist therapy against gastric cancer. Cancer Cell. 2014 Feb 10;25(2):166-80.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---