
Product Data Sheet

Product Name: Thioredoxin reductase peptide

Cat. No.: GC37782

Chemical Properties

Cas. No. 950890-23-0

Formula $C_{66}H_{106}N_{18}O_{18}S_2$ M.Wt 1503.79

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Thioredoxin reductase peptide corresponds to residues 53-67 in thioredoxin reductase (TrxR), used in thioredoxin reductase research. Thioredoxin reductase acts as a reductant of disulfide-containing proteins and plays crucial role in cellular antioxidant defense[1][2].

[1]. Jan YH, et al. Cross-linking of thioredoxin reductase by the sulfur mustard analogue mechlorethamine(methylbis(2-chloroethyl)amine) in human lung epithelial cells and rat lung: selective inhibition of disulfide reduction but not redox cycling. Chem Res Toxicol. 2014 Jan 21;27(1):61-75. [2]. Valette O, et al. Biochemical Function, Molecular Structure and Evolution of an Atypical Thioredoxin Reductase from *Desulfovibrio vulgaris*. Front Microbiol. 2017 Sep 29;8:1855.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
