
Product Data Sheet

Product Name: TLQP-30
Cat. No.: GC37800

Chemical Properties

Cas. No. 922704-13-0

Formula $C_{150}H_{232}N_{54}O_{40}$

M.Wt 3431.78

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

TLQP-30 is a VGF peptide[1].

[1]. D'Amato F, et al. Differential distribution of VGF-derived peptides in the adrenal medulla and evidence for their selective modulation. J Endocrinol. 2008 May;197(2):359-69.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA