
Product Data Sheet

Product Name: TNF- α (10-36), human (TFA)

Cat. No.: GC37804

Chemical Properties

Cas. No.

Formula C₁₃₃H₂₁₂N₄₃O₄₀ M.Wt 3110.36

Solubility Soluble in DMSO Storage -20°C, protect from light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**TNF- α (10-36), human (TFA) is a peptide of human TNF- α .

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA