

---

## Product Data Sheet

---

Product Name: NHPI-PEG4-C2-Pfp ester

Cat. No.: GC38018

### Chemical Properties

Cas. No. 2101206-46-4

SMILES O=C(OC1=C(F)C(F)=C(F)C(F)=C1F)CCOCCOCCOCCOCCON(C(C2=C3C=CC=C2)=O)C3=O

Formula C<sub>25</sub>H<sub>24</sub>F<sub>5</sub>NO<sub>9</sub>

M.Wt

577.45

Solubility Soluble in DMSO

Storage

Store at -20°C

General For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the tips ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Evaluation sample solution : ship with blue ice All other available size: ship with RT , or Condition blue ice upon request.

Structure

### Background

NHPI-PEG4-C2-Pfp ester is used as a linker for antibody-drug conjugates (ADC).

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA