

---

## Product Data Sheet

---

Product Name: Ald-Ph-amido-PEG3-C2-Pfp ester

Cat. No.: GC38072

### Chemical Properties

Cas. No. 2101206-21-5

SMILES O=C(OC1=C(F)C(F)=C(F)C(F)=C1F)CCOCCOCCOCCNC(C2=CC=C(C=O)C=C2)=O

Formula  $C_{23}H_{22}F_5NO_7$  M.Wt 519.42

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

Ald-Ph-amido-PEG3-C2-Pfp ester is a noncleavable ADC linker, which belongs to a polyethylene glycol (PEG) linker.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA