
Product Data Sheet

Product Name: Ald-Ph-amido-PEG3-C1-Boc

Cat. No.: GC38079

Chemical Properties

Cas. No. 1007215-94-2

SMILES O=C(OC(C)(C)C)COCCOCCOCCNC(C1=CC=C(C=O)C=C1)=OFormula C₂₀H₂₉NO₇ M.Wt 395.45

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Ald-Ph-amido-PEG3-C1-Boc is an ADC linker, which belongs to a polyethylene glycol (PEG) linker.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
