
Product Data Sheet

Product Name: PDB-Pfp
Cat. No.: GC38442

Chemical Properties

Cas. No. 2088570-81-2

SMILES O=C(OC1=C(F)C(F)=C(F)C(F)=C1F)CCCSSC2=NC=CC=C2

Formula $C_{15}H_{10}F_5NO_2S_2$ M.Wt 395.37

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

PDB-Pfp is a reducible ADC linker used in the synthesis of antibody-drug conjugates (ADCs).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA