
Product Data Sheet

Product Name: N3-PEG3-C2-NHS ester

Cat. No.: GC38455

Chemical Properties

Cas. No. 1245718-89-1

SMILES O=C(ON1C(CCC1=O)=O)CCOCCOCCOCCN=[N+]=[N-]Formula $C_{13}H_{20}N_4O_7$ M.Wt 344.32

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

N3-PEG3-C2-NHS ester is a noncleavable 3-unit PEG linker used in the synthesis of antibody-drug conjugates (ADCs).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
