

---

**Product Data Sheet**

---

Product Name: Ald-Ph-amido-PEG1-C2-NHS ester

Cat. No.: GC38476

**Chemical Properties**

Cas. No. 2101206-80-6

SMILES O=C(ON1C(CCC1=O)=O)CCOCCNC(C2=CC=C(C=O)C=C2)=OFormula  $C_{17}H_{18}N_2O_7$  M.Wt 362.33

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Ald-Ph-amido-PEG1-C2-NHS ester is a noncleavable 1-unit PEG linker for antibody-drug-conjugation (ADC).

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---