

---

**Product Data Sheet**

---

Product Name: 11-Aminoundecanoic acid

Cat. No.: GC38727

**Chemical Properties**

Cas. No. 2432-99-7

SMILES NCCCCCCCCCCC(O)=O

Formula C<sub>11</sub>H<sub>23</sub>NO<sub>2</sub>

M.Wt 201.31

Solubility DMSO: 3.33 mg/mL (16.54 mM; ultrasonic and adjust pH to 2 with HCl); Water: &lt; 0.1 mg/mL (insoluble)

Store  
Storage at -  
20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

11-Aminoundecanoic acid is a PROTAC linker, which refers to the alkyl chain composition. 11-Aminoundecanoic acid can be used in the synthesis of a series of PROTACs.

[1]. Pillow TH, et al. Antibody Conjugation of a Chimeric BET Degradable Enables in vivo Activity. ChemMedChem. 2019 Oct 31.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA