
Product Data Sheet

Product Name: Isoginsenoside Rh3

Cat. No.: GC39102

Chemical Properties

Cas. No. 166040-90-0

SMILES C[C@]12[C@@](C[C@@H](O)[C@@]3([H])[C@]2(CC[C@@H]3/C(C)=C/C/C=C(C)/C)C)([H])[C@@]4([C@@](C(C)([C@@H](O[C@]5([H])O[C@@H]([C@@H](O)[C@H](O)[C@H]5O)CO)CC4)C)([H])CC1)C

Formula C₃₆H₆₀O₇ M.Wt 604.86

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Isoginsenoside-Rh3 is a triterpenoid saponin from the fruits of Panax ginseng C. A. Mey[1].

[1]. Wang JY, et al. Isoginsenoside-Rh3, a new triterpenoid saponin from the fruits of Panax ginseng C. A. Mey. J Asian Nat Prod Res. 2004 Dec;6(4):289-93.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
