
Product Data Sheet

Product Name: Phosphorylcholine chloride calcium salt tetrahydrate

Cat. No.: GC39173

Chemical Properties

Cas. No. 72556-74-2

SMILES C[N+](C)(CCOP([O-])([O-])=O)C.[Ca+2].[Cl-].O.O.O.OFormula $C_5H_{21}CaClNO_8P$ M.Wt 329.73

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Phosphocholine (Phosphorylcholine) chloride calcium salt tetrahydrate is an active endogenous metabolite.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA