
Product Data Sheet

Product Name: Tos-PEG4-THP

Cat. No.: GC39207

Chemical Properties

Cas. No. 86259-89-4

SMILES O=S(C1=CC=C(C)C=C1)(OCCOCCOCCOCCOC2CCCCO2)=OFormula $C_{20}H_{32}O_8S$ M.Wt 432.53

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Tos-PEG4-THP is a PEG-based PROTAC linker can be used in the synthesis of PROTAC K-Ras Degradar-1 [1].

[1]. Andrew P. Crew, et al. Modulators of proteolysis and associated methods of use. WO2019195609A2.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
