
Product Data Sheet

Product Name: Amino-PEG4-C2-amine

Cat. No.: GC39216

Chemical Properties

Cas. No. 68960-97-4

SMILES NCCOCCOCCOCCOCCN

Formula C₁₀H₂₄N₂O₄ M.Wt 236.31

Solubility Soluble in DMSO Storage 4°C, protect from light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Amino-PEG4-C2-amine is a PEG-based (4 units) PROTAC linker can be used in the synthesis of PROTACs.

[1]. Lepage ML, et al. Design, synthesis and photochemical properties of the first examples of iminosugar clusters based on fluorescent cores. Beilstein J Org Chem. 2015 May 6;11:659-67.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
