
Product Data Sheet

Product Name: THP-PEG4-Pyrrolidine(N-Boc)-CH₂OH

Cat. No.: GC39228

Chemical Properties

Cas. No. 2378261-80-2

SMILES O=C(N1[C@H](CO)C[C@@H](OCCOCCOCCOCCOC2CCCCO2)C1)OC(C)(C)CFormula C₂₃H₄₃NO₉ M.Wt 477.59

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

THP-PEG4-Pyrrolidine(N-Boc)-CH₂OH is a PEG-based PROTAC linker can be used in the synthesis of PROTAC K-Ras Degradar-1 [1].

[1]. CREW, Andrew P., et al. MODULATORS OF PROTEOLYSIS AND ASSOCIATED METHODS OF USE. WO2019195609A2.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
