

---

**Product Data Sheet**

---

Product Name: Sulfo-DMAC-SPP

Cat. No.: GC39446

**Chemical Properties**

Cas. No. 663599-11-9

SMILES CC(SSC1=C(C(N(C)C)=O)C=NC=C1)CCC(ON2C(C(S(=O)(O)=O)CC2=O)=O)=OFormula C<sub>17</sub>H<sub>21</sub>N<sub>3</sub>O<sub>8</sub>S<sub>3</sub> M.Wt 491.56

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Sulfo-DMAC-SPP is a cleavable ADC linker used in the synthesis of antibody-drug conjugates (ADCs)[1].

[1]. Beck A, et al. Strategies and challenges for the next generation of antibody-drug conjugates. Nat Rev Drug Discov. 2017 May;16(5):315-337.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---