
Product Data Sheet

Product Name: Heptaethylene glycol

Cat. No.: GC39597

Chemical Properties

Cas. No. 5617-32-3

SMILES OCCOCCOCCOCCOCCOCCOCCO

Formula C₁₄H₃₀O₈ M.Wt 326.38

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Heptaethylene glycol (HO-PEG7-OH) is a PEG-based PROTAC linker that is applicable to the synthesis of PROTACs.

[1] Sainan An, Liwu Fu. EBioMedicine. 2018 Oct;36:553-562.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
