

---

## Product Data Sheet

---

Product Name: PBS, 100 Tablets

Cat. No.: GC50205

### Chemical Properties

Cas. No.

Formula

M.Wt

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

Phosphate buffered saline tablets. Dissolve 1 tablet into 100ml of water. This will produce 100ml of PBS solution containing: Sodium chloride (137mM), Potassium chloride (2.7mM), Phosphate buffer (11.9mM).

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA