

---

**Product Data Sheet**

---

Product Name: TP 472N  
Cat. No.: GC50312

**Chemical Properties**

Cas. No. 2080306-24-5

SMILES O=C(C)C1=CC(C2=CC=CC=C2OCC3CC3)=C4N=CC=CN41

Formula  $C_{19}H_{18}N_2O_2$  M.Wt 306.36

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

**Background**

Negative control for TP 472. Inactive against other bromodomains (>20 μM against BRD9). Active Analog also available.

0

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA