
Product Data Sheet

Product Name: H3K27(Me3) (15-34)

Cat. No.: GC50382

Chemical Properties

Cas. No.

Formula C₈₇H₁₅₇N₃₀O₂₅ M.Wt 2023.38

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Histone peptide. Substrate for methyltransferase/demethylase enzymes. Peptide derived from human histone isotype 3.1.

Chadwick and Willard et al (2004) Multiple spatially distinct types of facultative heterochromatin on the human inactive X chromosome. Proc.Natl.Acad.Sci.U.S.A. 101 17450 PMID:15574503

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
