
Product Data Sheet

Product Name: PEN (human)

Cat. No.: GC50416

Chemical Properties

Cas. No. 597578-70-6

Formula $C_{97}H_{159}N_{27}O_{32}$

M.Wt 2215.49

Solubility H_2O : 1.96 mg/mL (0.88 mM); ultrasonic and adjust pH to 3 with HCl

Storage Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Endogenous peptide GPR83 agonist. ProSAAS-derived neuropeptide.

Gomes et al (2016) Identification of GPR83 as the receptor for the neuroendocrine peptide PEN. Sci.Signal 9 ra43 PMID:27117253 |Mack et al (2019) Neuropeptide PEN and its receptor GPR83: distribution, signaling, and regulation. ACS Chem Neurosci. 10 1884 PMID:30726666

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
