
Product Data Sheet

Product Name: VH 032 amide-PEG3-acid

Cat. No.: GC50543

Chemical Properties

Cas. No. 2140807-42-5

SMILES CC1=C(C2=CC=C(C=C2)CNC([C@@H]3C[C@H](CN3C([C@H](C(C)(C)C)NC(CCOCCOCCOCCC(O)=O)=O)=O)O)=O)SC=N1

Formula	C ₃₂ H ₄₆ N ₄ O ₉ S	M.Wt	662.79
---------	---	------	--------

Solubility	Soluble in DMSO	Storage	Store at -20°C
------------	-----------------	---------	----------------

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Functionalized von-Hippel-Lindau (VHL) protein ligand for PROTAC® research and development; incorporates an E3 ligase ligand plus a PEG linker ready for conjugation to a target protein ligand. Part of a range of functionalized tool molecules for PROTAC R&D. This product has been recently renamed. The previous name for this product was VH 032 - linker 4 PROTAC® is a registered trademark of Arvinas Operations, Inc., and is used under license.

0

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
