
Product Data Sheet

Product Name: Amino-PEG3-C2-acid

Cat. No.: GC60046

Chemical Properties

Cas. No. 784105-33-5

SMILES O=C(O)CCOCCOCCOCCNFormula $C_9H_{19}NO_5$

M.Wt

221.25

Solubility

Storage

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Amino-PEG3-C2-acid is a cleavable PEG ADC linker used in the synthesis of antibody-drug conjugates (ADCs)[1].

[1]. Kevin HAMBLETT, et al. Anti-her2 biparatopic antibody-drug conjugates and methods of use. WO2019173911A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
