
Product Data Sheet

Product Name: Azido-PEG3-SS-NHS

Cat. No.: GC60068

Chemical Properties

Cas. No.

SMILES O=C(ON1C(CCC1=O)=O)CCSSCCOCCOCCOCCN=[N+]=[N-]Formula C15H24N4O7S2 M.Wt 436.5

Solubility Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Azido-PEG3-SS-NHS is a cleavable 3 unit PEG ADC linker used in the synthesis of antibody-drug conjugates (ADCs)[1].

ADCs are comprised of an antibody to which is attached an ADC cytotoxin through an ADC linker[1].

[1]. Beck A, et al. Strategies and challenges for the next generation of antibody-drug conjugates. Nat Rev Drug Discov. 2017 May;16(5):315-337.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA