
Product Data Sheet

Product Name: Propargyl-PEG3-alcohol

Cat. No.: GC60304

Chemical Properties

Cas. No. 208827-90-1

SMILES C#CCOCCOCCOCCOFormula $C_9H_{16}O_4$

M.Wt

188.22

Solubility

Storage

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Propargyl-PEG3-alcohol is a PEG-based PROTAC linker can be used in the synthesis of PROTACs[1].

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins.

[1]. Andrew P. Crew, et al. Brm targeting compounds and associated methods of use. US20190300521A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
