
Product Data Sheet

Product Name: 2-Acetyl-3-ethylpyrazine

Cat. No.: GC60469

Chemical Properties

Cas. No. 32974-92-8

SMILES CC(C1=NC=CN=C1CC)=OFormula C₈H₁₀N₂O

M.Wt 150.18

Solubility DMSO : 100 mg/mL (665.87 mM; Need ultrasonic)

Storage 4°C, stored under nitroge

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

2-Acetyl-3-ethylpyrazine is an endogenous metabolite.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA