

---

**Product Data Sheet**

---

Product Name: 3,5-Dihydroxyacetophenone

Cat. No.: GC60497

**Chemical Properties**

Cas. No. 51863-60-6

SMILES CC(C1=CC(O)=CC(O)=C1)=OFormula C<sub>8</sub>H<sub>8</sub>O<sub>3</sub>

M.Wt 152.15

Solubility DMSO : 50 mg/mL (328.62 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

3,5-Dihydroxyacetophenone is an endogenous metabolite.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA