
Product Data Sheet

Product Name: Mid-chain triglyceride of caprylic/caprolic acid

Cat. No.: GC61063

Chemical Properties

Cas. No.

SMILES [Mid-chain triglyceride of caprylic/caprolic acid]

Formula M.Wt

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Mid-chain triglyceride of caprylic/caprolic acid can be used as a co-solvent.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA