

---

**Product Data Sheet**

---

Product Name: Tiopronin 13C D3

Cat. No.: GC61334

**Chemical Properties**

Cas. No. 1189695-13-3

SMILES O=C(O)[13C]NC(C(S)C([2H])([2H])[2H])=OFormula  $C_4\ 13CH_4D_3NO_3S$  M.Wt 165.19

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Tiopronin 13C D3 is deuterium labeled Tiopronin.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA