
Product Data Sheet

Product Name: DBCO-PEG4-MMAF

Cat. No.: GC62261

Chemical Properties

Cas. No. 2360411-65-8

Formula C₆₉H₉₉N₇O₁₅

M.Wt 1266.56

Solubility DMSO : 200 mg/mL (157.91 mM;
Need ultrasonic)Storage Store at -20°C, unstable in
solution, ready to use.

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

DBCO-PEG4-MMAF is a drug-linker conjugate for ADC with potent antitumor activity by using the tubulin polymerization inhibitor, MMAF, linked via the cleavable linker DBCO-PEG4.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA