
Product Data Sheet

Product Name: Fmoc-D-Val-D-Cit-PAB

Cat. No.: GC62373

Chemical Properties

Cas. No. 1350456-69-7

Formula $C_{33}H_{39}N_5O_6$ M.Wt 601.69Solubility Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Fmoc-D-Val-D-Cit-PAB is a cleavable linker for antibody-drug-conjugation (ADC)[1].

[1]. Peter D. Senter, et al. Drug conjugates and their use for treating cancer, an autoimmune disease or an infectious disease. WO2004010957A2.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
